

Bioconcept
Bone Level (Tapered) Implant System

Optimize crestal bone preservation and soft tissue stability
Apically tapered implant body design allows under-preparation and
supports a high primary stability in soft bone



Benefit more people



Protective Cap					
	Diameter		Φ3.5mm	Diameter	Φ4.6mm
	H	5.0mm	123010	H	5.1mm 123040
		6.5mm	123020		6.6mm 123050
		8.0mm	123030		8.1mm 123060

* Recommended screw tightening torque: 15Ncm

Protective Cap			
	Diameter		Φ4.6mm
	H	5.1mm	123040
		6.6mm	123050
		8.1mm	123060

* Recommended screw tightening torque: 15Ncm

Manual Impression				Digital Impression			
Closed Tray		Open Tray		Scanbody		Analog	
	Φ3.5		Φ3.5		Φ3.5		Φ3.5
	Crown		Crown		Crown		0°
	Bridge		Bridge		Bar/Bridge		162750
	Crown		Crown		Crown		0°
	162730		162710		162150		17°/30°
	Bridge		Bridge		162710N		162760
	162730N				162150N		

Manual Impression				Digital Impression			
Closed Tray		Open Tray		Scanbody		Analog	
	Φ4.6		Φ4.6		Φ4.6		Φ4.6
	Crown		Crown		Crown		
	Bridge		Bridge		Bridge		
	162730		162710		162150		
	162730N		162710N		162150N		

Screw-retained

Screw-retained Abutment (including basal screw)							
	↖		0°		17°		30°
	Φ		3.5mm	4.6mm	4.6mm	4.6mm	
	GH	1mm	131410	131440	-	-	
		2.5mm	131420	131450	Type A	131470	Type A 131530
					Type B	131480	Type B 131540
		4mm	131430	131460	Type A	131490	Type A 131550
					Type B	131500	Type B 131560
	5.5mm	-	-	-	Type A	131510	Type A 131570
					Type B	131520	Type B 131580

* Recommended screw tightening torque: 35Ncm

Screw-retained Abutment (including basal screw)							
	↖		0°		17°		30°
	Φ		4.6mm	4.6mm	4.6mm	4.6mm	
	GH	1mm	132410	-	-	-	
		2.5mm	132420	Type A	132440	Type A	132500
				Type B	132450	Type B	132510
		4mm	132430	Type A	132460	Type A	132520
				Type B	132470	Type B	132530
	5.5mm	-	-	Type A	132480	Type A	132540
				Type B	132490	Type B	132550

* Recommended screw tightening torque: 35Ncm

Titanium Temporary Coping				TiBase Abutment for Bridge/Bar Cylindrical Coping			
	Φ	3.5mm	4.6mm		Φ3.5mm	Φ4.6mm	
	Crown	171410	171420		811170N	811180N	
	Bridge	171410N	171420N		including screw		

* Recommended screw tightening torque: 15Ncm

* Recommended screw tightening torque: 15Ncm

Titanium Temporary Coping			TiBase Abutment for Bridge/Bar Cylindrical Coping		
	Φ	4.6mm		Φ4.6mm	
	Crown	171420		811180N	
	Bridge	171420N		including screw	

* Recommended screw tightening torque: 15Ncm

* Recommended screw tightening torque: 15Ncm

Occlusal Screw for Coping, Screw-retained Abt		Lab Processing Screw for Multi Abt		Polishing Protector for Multi Abt	
	141010		L20 164010 L10 164020		Φ3.5 157090 Φ4.6 157100

* Recommended screw tightening torque: 15Ncm

Occlusal Screw for Coping, Screw-retained Abt		Lab Processing Screw for Multi Abt		Polishing Protector for Multi Abt	
	141010		L20 164010 L10 164020		Φ4.6 157100

* Recommended screw tightening torque: 15Ncm

Abutment Level Impression		
Equator		

Impression Coping		
	368010	
Analog		
	368020	





Ball-retained

Equator Abutment			
	GH	NC	RC
	1		134210
	2	134020	134220
	3	134030	134230
	4	134040	134240
	5	134050	134250
	6	134060	134260

* Recommended screw tightening torque: 35Ncm

SmartBox Male Processing Package 378100					
Processing Male Protective Disk		Equator Replacement Male			
		2. 7kg	1. 8kg	1. 2kg	0. 6kg
378070v20		378020v4	378030v4	378040v4	378050v4

Equator Square Screwdriver for ratchet	
	055130
Equator Insertion-Extraction Tool	
	358020

Bioconcept		Φ3.3mm		Φ4.1mm		Φ4.8mm				NC		RC	
Bioconcept Bone Level Implant System Tel: +86-519-67895691 E-mail: intlsales@bioconcept.cn https://www.bioconcept.cn Add: No.1019, Liaohe Road, Xinbei District, Changzhou, Jiangsu, 213125, China		B-L Implant NC Φ3.3mm  Φ3.3mm	B-L-T Implant NC Φ3.3mm 	B-L Implant RC Φ4.1mm  Φ4.1mm	B-L-T Implant RC Φ4.1mm 	B-L Implant RC Φ4.8mm  Φ4.8mm	B-L-T Implant RC Φ4.8mm 			Implant Level Impression		Implant Level Impression	
		Manual Impression		Digital Impression		Manual Impression		Digital Impression					
		NC Impression Post for Closed Tray 161090 NC Impression Post Open Tray 161100 / 161110   		NC Scanbody 861010 		RC Impression Post for Closed Tray 162090 RC Impression Post Open Tray 162100 / 162110   		RC Scanbody 861020 					
		NC Fixture Analog 161010 		NC Digital Analog 161800 (Φ3.6mm) 		RC Fixture Analog 162010 		RC Digital Analog 162820 (Φ4.8mm) 					